



Balanzas de precision



professional level



PS/Y

- max: 250 g / 0,001 g
- max: 7800 g / 0,01 g
- graphic display 5,7" with touch panel
- automatic internal calibration



WLY

- max: 3/6 kg / 0,05/0,1 g
- max: 60/120 kg / 1/2 g
- graphic display 5,7" with touch panel



APP/Y

- max: 35 kg / 0,1 g
- graphic display 5,7" with touch panel
- automatic internal calibration

standard level



PS/X

- max: 200 g / 0,001 g
- max: 7800 g / 0,01 g
- graphic display

basic level



PS /C1

- max: 110 g / 0,001 g
- max: 6000 g / 0,01 g
- LCD display



PS /C2

- max: 110 g / 0,001 g
- max: 6000 g / 0,01 g
- LCD display
- automatic internal calibration



APP/C

- max: 6/35 kg
- LCD display
- electromagnetic measuring mechanism



WLC

- max: 0,6 kg / 0,01 g
- max: 120 kg / 2 g
- LCD display



WLC/C2

- max: 0,6 kg / 0,01 g
- max: 6 kg / 0,1 g
- LCD display



WTB

- max: 3000 g / 0,01 g
- max: 200 g / 0,001 g
- LCD display

MICROSERCON,SLU; CIF: B-83691113

Urb.Rio Cofio C/Los Pinos 217 28294 – Robledo de Chavela (Madrid)

Tlf 918998029 / Fax 918998682 Dpto.Comercial 677590684 Dpto.Tecnico 678558460

[E mail: microsercon@telefonica.net](mailto:microsercon@telefonica.net) www.microsercon.com