

**CAMARAS INFINITY****PVP €****INFINITY LiteB CMOS USB 2,0 Camera**

<b>INFINITY LiteB</b>	Sensor 1/2,5"CMOS / Pixels 1.5 Mp / Frame Rates: 15 Resolución: 1400 x 1080 / Square pixels 4,2 µm	<b>1.100,00 €</b>
-----------------------	---	-------------------

**INFINITY 1 CMOS USB 2,0 Camera**

<b>INFINITY 1-2 CB</b>	Sensor 1/2"CMOS / Pixels 2,0 Mp / Frame Rates: 15 Resolución: 1600 x 1200 / Square pixels 4,2 µm	<b>1.600,00 €</b>
<b>INFINITY 1-3 C</b>	Sensor 1/2"CMOS / Pixels 3,1 Mp / Frame Rates: 12 Resolución: 2048 x 1536 / Square pixels 3,2 µm	<b>2.400,00 €</b>
<b>INFINITY 1-5 C</b>	Sensor 1/2,5"CMOS / Pixels 5,0 Mp / Frame Rates: 5 Resolución: 2592 x 1944 / Square pixels 2,2 µm	<b>3.600,00 €</b>

**INFINITY 2 CCD USB 2,0 Camera**

<b>INFINITY 2-1 C/M</b>	Sensor 1/2"CCD / Pixels 1,4 Mp / Frame Rates: 15 Resolución: 1392 x 1040 / Square pixels 4,65 µm	<b>2.200,00 €</b>
<b>INFINITY 2-2 C/M</b>	Sensor 1/2"CCD / Pixels 2,0 Mp / Frame Rates: 12 Resolución: 1616 x 1216 / Square pixels 4,4 µm	<b>3.300,00 €</b>
<b>INFINITY 2-3 C/M</b>	Sensor 1/8"CCD / Pixels 3,3 Mp / Frame Rates: 5 Resolución: 2080 x 1536 / Square pixels 3,45 µm	<b>3.600,00 €</b>
<b>INFINITY 2-5 C/M</b>	Sensor 1/8"CCD / Pixels 5,0 Mp / Frame Rates: 8 Resolución: 2080 x 1536 / Square pixels 3,45 µm	<b>4.500,00 €</b>

**INFINITY 3 CCD USB 2,0 Camera**

<b>INFINITY 3-1 C/M</b>	Sensor 2/3"Cooled CCD / Pixels 1,4 Mp / Frame Rates: 15 Resolución: 1392 x 1040 / Square pixels 6,45 µm	<b>7.100,00 €</b>
<b>INFINITY 3-1U C/M</b>	Sensor 2/3" CCD / Pixels 1,4 Mp / Frame Rates: 15 Resolución: 1392 x 1040 / Square pixels 6,45 µm	<b>5.200,00 €</b>

**INFINITY 4 CCD USB 2,0 Camera**

<b>INFINITY 4-11 C/M</b>	Sensor 35mm Format CCD / Pixels 10,7 Mp / Frame Rates: 3 Resolución: 4008 x 2670 / Square pixels 9,0 µm	<b>8.200,00 €</b>
--------------------------	--	-------------------

**INFINITY X CCD USB 2,0 Camera**

<b>INFINITY X-21 M</b>	Sensor Pixel Shifting 1/2" CCD / Pixels 21 Mp / Frame Rates: 15 Resolución: up to 5120 x 4096 / Square pixels 4,4 µm	<b>4.100,00 €</b>
<b>INFINITY X-32 C/M</b>	Sensor Pixel Shifting 1/1,8" CCD / Pixels 32 Mp / Frame Rates: 12 Resolución: 6484 x 4864 / Square pixels 4,4 µm	<b>6.300,00 €</b>

**ADVANCED FEATURES SOFTWARE MODULE****850,00 €****MICROSERCON,SLU**; CIF: B-83691113

Urb.Rio Cofio C/Los Pinos 217 28294 – Robledo de Chavela (Madrid)

Tif 918998029 / Fax 918998682 Dpto.Comercial 677590684 Dpto.Tecnico 678558460

[E mail: microsercon@telefonica.net](mailto:microsercon@telefonica.net) [www.microsercon.com](http://www.microsercon.com)